



Buzsaki256 / A8x32(linear)


p.1

127	98	81	67	206	224	239	241
115	103	88	76	192	213	230	253
123	100	83	71	194	222	237	250
124	105	90	66	195	211	228	249
119	102	85	75	198	220	235	246
120	107	92	70	199	209	226	245
114	104	87	78	202	218	233	242
116	109	94	74	203	207	223	255
177	188	184	180	62	58	54	51
191	187	183	179	49	61	57	53
168	174	172	170	42	44	46	48
167	165	163	161	32	33	35	37
151	157	155	153	25	27	29	31
152	150	148	146	16	18	20	22
142	132	136	140	0	4	8	12
128	131	135	139	14	5	9	13
141	129	133	137	3	7	11	1
130	134	138	144	15	2	6	10
143	149	147	145	17	19	21	23
160	158	156	154	24	26	28	30
159	166	164	162	34	36	38	40
175	173	171	169	39	41	43	45
176	189	185	181	63	59	55	47
190	186	182	178	50	60	56	52
117	101	86	77	204	215	232	254
113	112	95	79	201	210	225	243
121	99	84	73	200	217	234	244
118	110	93	72	197	212	227	247
125	97	82	69	196	219	236	248
122	108	91	68	193	214	229	251
111	96	80	65	208	221	238	252
126	106	89	64	205	216	231	240



A4x64-Poly2-5mm-23s-250-177

p.2





A4x64-Poly2-10mm-23s-200-177


p.3

142		131	136		139	0		5	8		13
151		150	155		146	25		18	29		22
168		165	172		161	42		33	46		37
177		187	184		179	62		61	54		53
114		109	87		74	202		207	233		255
119		107	85		70	198		209	235		245
123		105	83		66	194		211	237		249
127		103	81		76	206		213	239		253
115		98	88		67	192		224	230		241
124		100	90		71	195		222	228		250
120		102	92		75	199		220	226		246
116		104	94		78	203		218	223		242
191		188	183		180	49		58	57		51
167		174	163		170	32		44	35		48
152		157	148		153	16		27	20		31
128		132	135		140	14		4	9		12
130		129	138		137	15		7	6		1
160		149	156		145	24		19	28		23
175		166	171		162	39		36	43		40
190		189	182		181	50		59	56		47
113		101	95		77	201		215	225		254
118		99	93		73	197		217	227		244
122		97	91		69	193		219	229		248
126		96	89		65	205		221	231		252
141		134	133		144	3		2	11		10
143		158	147		154	17		26	21		30
159		173	164		169	34		41	38		45
176		186	185		178	63		60	55		52
117		112	86		79	204		210	232		243
121		110	84		72	200		212	234		247
125		108	82		68	196		214	236		251
111		106	80		64	208		216	238		240



A8x32-Edge-5mm-25-200-177


p.4





A8x32-Poly2-5mm-20s-150-160


p.5





A8x32-Poly2-5mm-20s-150-160


p.6






A8x32-Poly2-5mm-20s-150-160

p.7



A8x32-Poly2-1.5mm-25s-325-177


p.8





A8x32-Poly2-1.5mm-25s-325-177

p.9





A8x32-Poly2-1.5mm-25s-325-177

p.10

